



QUIDEL®

QuickVue® At-Home OTC COVID-19 Test

USER INSTRUCTIONS


For Emergency Use
Authorization (EUA) only.
In vitro diagnostic use only.



Scan QR Code with phone camera to learn
more or visit quickvueathome.com.

Wash Your Hands


Before you start testing,
wash your hands or use
hand sanitizer.
Make sure your hands are
dry before starting.



STEP 1

Check Your Test Kit

Locate the kit components:
It is recommended gloves (not provided)
also be used during testing.



Tubes



Tube Holder



Swabs



Watch or Timer
(not included)




Test Strips

STEP 2

Place Tube in the Tube Holder

Remove cap from one
TUBE and place it in the
TUBE HOLDER.


NOTE: Use of gloves
is recommended.



STEP 3

Swab the Nostrils

A Remove the **SWAB** from
its wrapper, being careful
not to touch the
SWAB head.




B Gently insert the **SWAB** ½ to ¾ of an inch into the nostril, depending on the
size of the person's nose. Firmly rub the **SWAB** in a circular motion around
the inside wall of **EACH NOSTRIL** at least **4** times.

Be sure to rub BOTH nostrils with the SAME SWAB.

NOTE: If you are swabbing others, please wear a face mask.
With children, you may not need to insert the swab as far into
the nostril. For very young children,
you may need another person
to steady the child's head
while swabbing.


NOTE: Failure to swab
properly may cause
false negative results.



STEP 4


Place Swab in the Tube

A Immediately place
the **SWAB** into the
liquid inside the
TUBE, and ensure
it is touching
the bottom.
Stir 3-4 times.



B Leave the swab in the
solution for **ONE MINUTE**.


NOTE: If the swab is in
the solution for more
than 10-minutes it
should not be used.




STEP 5

Remove Swab from the Tube

A After **ONE MINUTE**, remove
the swab from the **TUBE** by
rubbing the swab
head against the
inside wall of the
tube to squeeze
out as much
liquid as possible.



B Dispose of swab
in the trash.




STEP 6

Open the Test Strip

Open the **TEST STRIP** pouch
carefully at the slit and hold
the **TEST STRIP** as indicated.




▼ Tear here



STEP 7

Place Test Strip in the Tube

Place the **TEST STRIP** into the **TUBE**
with the arrows pointing down.
Leave the strip in the **TUBE**
for a **FULL TEN MINUTES** -
do not handle or remove.



Leave strip
in at least
10 minutes!

Arrows
pointing
down

STEP 8

Remove Test Strip from the Tube

At **TEN MINUTES**, remove the **TEST STRIP** from
the **TUBE**. Next, place the **TEST STRIP** on the
outline in Steps 10-12. Ensure the **TEST STRIP**
is on a flat surface in good lighting.

NOTE: The test is intended to be read at
10 minutes. If the test is read before this
or is read more than 5 minutes after the
indicated read time, results may be
inaccurate (false negative, false
positive, or invalid) and the
test should be repeated.



STEP 9

Check Your Results

There are three type of
results possible.

1. Check for a Positive Result
2. Check for a Negative Result
3. Check for an Invalid Result

STEP 10

Check for a Positive COVID-19 Result


Place the **TEST STRIP** on the test strip outline below and compare
with test result examples shown.




Positive Result

A **POSITIVE** result must show
BOTH a **BLUE** line and a **PINK**
line near the **BLUE** line.

Look closely! Even a very faint,
pink Test Line and a blue
Control Line is a **POSITIVE**
result. The intensity of the
lines may vary.



Results shown at 2x.



Positive COVID-19 Result

A positive test result means that proteins from the virus
that causes COVID-19 were found in your sample and it
is very likely you have COVID-19 and it is important to be
under the care of your healthcare provider. It is also likely
that you may be placed in isolation to avoid spreading
the virus to others. There is a very small chance that this
test can give you a positive test result that is wrong (false
positive.) If you test positive with the QuickVue At-Home
OTC COVID-19 Test you should self-isolate and seek
follow-up care with your healthcare provider as
additional testing may be necessary. Your healthcare
provider will work with you to determine how best to
care for you based on your test result along with your
medical history, and your symptoms.



QUIDEL®

QuickVue® At-Home OTC COVID-19 Test

USER INSTRUCTIONS

For Emergency Use
Authorization (EUA) only.
In vitro diagnostic use only.

Continued from other side

STEP 11

Check for a Negative COVID-19 Result

Place the **TEST STRIP** on the test strip outline below and compare with test result examples shown.



Negative Result

A **NEGATIVE** result will show a **BLUE** line but **NO PINK** line.



Results shown at 2x.



Negative COVID-19 Result

A negative test result means that proteins from the virus that causes COVID-19 were not found in your sample.

It is possible for this test to give a negative result that is incorrect (false negative) in some people with COVID-19. This means that you could possibly still have COVID-19 even though the test is negative. If you test negative and continue to experience COVID-19 like symptoms of fever, cough and/or shortness of breath you should seek follow up care with your healthcare provider. Your healthcare provider will consider the test result together with all other aspects of your medical history (such as symptoms, possible exposures, and geographical location of places you have recently traveled) in deciding how to care for you. For example, your healthcare provider may suggest you need another test to determine if you have contracted the virus causing COVID-19.

It is important that you work with your healthcare provider to help you understand the next steps you should take.

STEP 12

Check for an Invalid COVID-19 Result

Place the **TEST STRIP** on the test strip outline below and compare with test result examples shown.



Invalid Result

If there is **NO LINE**, or if there is **ONLY a PINK** line, the test is **INVALID** and you should repeat the steps starting at the beginning using a new **TUBE**, **SWAB**, and **TEST STRIP**.




Results shown at 2x.



STEP 13

Dispose Used Test in the Trash

All used test components should be disposed of in your household waste.



Wash Your Hands

After completing all steps, wash hands or use hand sanitizer.



Please notify your healthcare provider of the results of your QuickVue At-Home OTC COVID-19 Test.

A second test should be obtained with at least 24 hours (and no more than 36 hours) between tests.

The QuickVue At-Home OTC COVID-19 Test is for FDA Emergency Use Authorization (EUA) Only For In Vitro Diagnostic (IVD) Use.

- This product has not been FDA cleared or approved but has been authorized by FDA under an EUA, not for any other viruses or pathogens.
- This product has been authorized only for the detection of proteins from SARS-CoV-2.
- The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of IVDs for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.
- For more information on EUAs go here: <https://www.fda.gov/emergency-preparedness-and-response-mcm-legal-regulatory-and-policy-framework/emergency-use-authorization>
- For the most up to date information on COVID-19, please visit: www.cdc.gov/COVID19
- For detailed instructions, please visit www.quickvueathome.com

Warnings, Precautions and Safety Information

- The QuickVue At-Home OTC COVID-19 Test is intended for the qualitative detection of the nucleocapsid protein antigen from SARS-CoV-2 from individuals with or without symptoms of COVID-19. The test is intended to be used twice over two to three days, with at least 24 hours and no more than 36 hours between tests.
- Read the written instructions fully before starting the test procedure
- To ensure correct results, you must follow the instructions
- Keep test kit and materials out of the reach of children and pets before and after use
- Wear safety mask or other face covering when collecting swabs from children or others
- Use of personal protection materials such as gloves are recommended
- Do not open the materials until ready for use. If the test strip is open for an hour or longer, invalid test results may occur.
- Improper swab collection may result in incorrectly negative (false negative) results
- The test is intended to be read at 10 minutes. If the test is read before this or is read more than 5 minutes after the indicated read time, results may be inaccurate and the test should be repeated.
- Do not use a test kit that is expired
- Do not touch the swab head when handling the swab
- Avoid exposure of your skin, eyes, nose, or mouth to the solution in the tube.

Hazardous Ingredients for Liquid Reagent

Chemical Name/CAS	Hazms (GHS Code) for each ingredient	Concentration
Sodium Phosphate Monobasic Monohydrate/10049-21-5	Causes skin irritation (H315) Causes serious eye irritation (H319) May cause respiratory irritation (H335)	0.7%
Sodium Phosphate Dibasic Anhydrous/7558-79-4	Causes serious eye damage (H318) Causes serious eye irritation (H319)	0.7%
C12-14-Alkylidimethyl betaines/66455-29-6	Causes severe skin burns and eye damage (H314) Causes serious eye damage (H318) Causes skin irritation (H315) Causes serious eye irritation (H319)	0.03%
ProClin® 300	Harmful if swallowed (H302) Harmful if inhaled (H332) May cause an allergic skin reaction (H317)	0.03%
ED1A Tetrasodium Salt/64-02-8	Harmful if swallowed (H302) Causes serious eye damage (H318) Causes serious eye irritation (H319) Harmful if inhaled (H332) May cause respiratory irritation (H335) May cause damage to organs (H371), single exposure	0.2%

- The solution in the tube contains hazardous ingredients (see table above). If the solution contacts the skin or eye, flush with plenty of water. If irritation persists, seek medical advice: <https://www.poison.org/contact-us-or-1-800-222-1222>

Intended Use

The QuickVue At-Home OTC COVID-19 Test is intended for the qualitative detection of the nucleocapsid protein antigen from SARS-CoV-2 from individuals with or without symptoms or other epidemiological reasons to suspect COVID-19 when tested twice over two or three days with at least 24 hours and no more than 36 hours between tests. This test is authorized for non-prescription home use with self-collected (unobserved) direct anterior nasal (NS) swab specimens from individuals aged 14 years and older or with adult-collected anterior NS samples from individuals aged 2 years or older.

The QuickVue At-Home OTC COVID-19 Test does not differentiate between SARS-CoV and SARS-CoV-2. Results are for the identification of SARS-CoV-2 nucleocapsid protein antigen. Antigen is generally detectable in anterior nasal specimens during the acute phase of infection. Positive results indicate the presence of viral antigens, but clinical correlation with past medical history and other diagnostic information is necessary to determine infection status. Positive results do not rule out bacterial infection or co-infection with other viruses and the agent detected may not be the definite cause of disease. Individuals who test positive with the QuickVue At-Home OTC COVID-19 Test should self-isolate and seek follow-up care with their physician or healthcare provider as additional testing may be necessary. Negative results should be treated as presumptive, do not rule out SARS-CoV-2 infection and should not be used as the sole basis for treatment or patient management decisions. In addition, negative results should be considered in the context of an individual's recent exposures, history and the presence of clinical signs and symptoms consistent with COVID-19, and confirmed with a molecular assay, if necessary, for patient management. For serial testing programs, additional confirmatory testing with a molecular test for negative results may be necessary, if there is a high likelihood of SARS-CoV-2 infection, such as in an individual with as a close contact with COVID-19 or with suspected exposure to COVID-19 or in communities with high prevalence of infection. Additional confirmatory testing with a molecular test for positive results may also be necessary, if there is a low likelihood of SARS-CoV-2 infection, such as in individuals without known exposures to SARS-CoV-2 or residing in communities with low prevalence of infection.

Individuals who test negative and continue to experience COVID-19 like symptoms of fever, cough and/or shortness of breath may still have SARS-CoV-2 infection and should seek follow up care with their physician or healthcare provider. Individuals should provide all results obtained with this product to their healthcare provider for public health reporting. All healthcare providers will report all test results they receive from individuals who use the authorized product to relevant public health authorities in accordance with local, state, and federal requirements using appropriate IONC and SNOMED codes, as defined by the Laboratory In Vitro Diagnostics (LIVD) Test Code Mapping for SARS-CoV-2 Tests provided by CDC.

The QuickVue At-Home OTC COVID-19 Test is authorized for non-prescription self-use and/or, as applicable for an adult lay user testing, another person aged 2 years or older in a non-laboratory setting. The QuickVue At-Home OTC COVID-19 Test is only for use under the Food and Drug Administration's Emergency Use Authorization.

What do I need to know about Results from Serial Testing?

If your first or second test is positive, then proteins from the virus that causes COVID-19 has been found in your sample and you likely have COVID-19. If you test positive with the QuickVue At-Home COVID-19 Test, you should self-isolate and seek follow-up care with your healthcare provider to determine the next steps you should take. You may need additional testing, depending on your personal health history and other factors. If your first test is negative, you should test again in 24-36 hours. If both your first and second tests are negative, you may not have COVID-19, however, you should follow-up with your healthcare provider if you are at high risk for COVID-19 infection. If you test negative and continue to experience COVID-19 like symptoms of fever, cough and/or shortness of breath you should seek follow up care with your healthcare provider. For example, your healthcare provider may suggest you need another test to determine if you have contracted the virus causing COVID-19. It is important that you work with your healthcare provider to help you understand the next steps you should take.

Frequently Asked Questions

Will this test hurt?

No, the nasal swab is not sharp and it should not hurt. Sometimes the swab can feel slightly uncomfortable. If you feel pain, please stop the test and seek advice from a healthcare provider.

What are the known and potential risks and benefits of this test?

Potential risks include:

- Possible discomfort during sample collection.
- Possible incorrect test results (see Results section).

Potential benefits include:

- The results, along with other information, can help your healthcare provider make informed recommendations about your care.
- The results of this test may help limit the spread of COVID-19 to your family and others in your community.

What is Serial Testing?

COVID-19 Serial testing is when one person tests themselves multiple times for COVID-19 on a routine basis, such as every day or every other day. By testing more frequently, you may detect COVID-19 more quickly and reduce spread of infection.

What is the difference between an antigen and molecular test?

An antigen test, such as the QuickVue At-Home OTC COVID-19 Test, detects proteins from the virus. Molecular tests detect genetic material from the virus. Antigen tests are very specific for the virus, but not as sensitive as molecular tests. This means that a positive result is highly accurate, but a negative result does not rule out infection. If your test result is negative, you should discuss with your healthcare provider on whether an additional test is necessary and if you should continue isolating at home.

How Accurate is this Test?

Based on the interim results of a clinical study where the QuickVue At-Home OTC COVID-19 Test was compared to an FDA authorized molecular SARS-CoV-2 test, QuickVue At-Home OTC COVID-19 Test correctly identified 83.5% of positive specimens and 99.2% of negative specimens. The performance of this test is still being studied in patients without signs and symptoms of respiratory infection and for serial screening. Performance may differ in these populations.